







DEVICE SPECS

ENDOPATH™ Endoscopic Hand Instruments

	 ENDOPATH™ Endoscopic Dissector	 ENDOPATH™ Endoscopic Blunt-Tip Dissector	 ENDOPATH™ Endoscopic Scissors	 ENDOPATH™ Endoscopic Grasper	 ENDOPATH™ Endoscopic Anvil Grasper	 ENDOPATH™ Endoscopic Babcock
Codes	5DCD	BTD05, BCD10	5DCS	5DSG	10AG	5BB, 10BB
Description	Dissector, curved, with monopolar cautery	Dissector, blunt tip/blunt tip (cherry)	Scissors, curved, with monopolar cautery	Grasper	Anvil Grasper	Babcock
Shaft Diameter	5mm	5mm or 10mm	5mm	5mm	10mm	5mm or 10mm
Shaft Rotation	360°	360°	360°	360°	360°	360°
Qty./Box	6	36 (3 per pouch, 12 pouches per sales unit)	6	6	6	6