

DEVICE SPECS

ENDOPATH™ XCEL Trocar Portfolio



ENDOPATH XCEL™ Trocars with OPTIVIEW™ Technology



ENDOPATH XCEL™ Trocars

	ENDOPATH XCEL™ Trocars with OPTIVIEW™ Technology	ENDOPATH XCEL™ Trocars
Codes	2B5ST, 2B5LT, 2B5XT, 2B12LT, 2B12XT, 2H12LP, 2D5ST, 2D5LT, 2D12LT, 2D12XT, 2CB5ST, 2CB5LT, 2CB12LT	B5ST, B5LT, B5XT, B8LT, B11LT, B11LP, B11LTH, B11LPH, B12SRT, B12LT, B12LP, B12LTH, B12LPH, B12XT, B15LT, D5ST, D5LT, D11LT, D12LT, D12XT, H12LP, CB5ST, CB5LT, CB11LT, CB12LT
Port Application	Camera Ports	Working Ports
Available Tip Configurations	Bladeless, Blunt Tip, Dilating Tip	Bladeless, Dilating Tip, Blunt Tip
Available Diameters	5mm, 12mm	5mm, 8mm, 11mm, 12mm, 15mm
Available Lengths	75mm, 100mm, 150mm	75mm, 100mm*, 150mm
Available Sleeve Types	Stability, Smooth	Stability, Smooth
Qty./Box	6	6