

DEVICE SPECS

# ENDOPOUCH™ Specimen Retrieval Bag



**ENDOPOUCH™**  
Specimen Retrieval Bag

<b>Codes</b>	POUCH
<b>Description</b>	Specimen retrieval bag
<b>Shaft Diameter</b>	10mm
<b>Shaft Length</b>	262mm
<b>Pouch Size</b>	102mm x 152mm (4in x 6in)
<b>Pouch Volume</b>	224mL
<b>Pouch Material</b>	Polyurethane Elastomer
<b>Deployment</b>	One-handed syringe motion
<b>Qty./Box</b>	6