

DEVICE SPECS

Ethicon Hernia Mesh Fixation Devices



ETHICON SECURESTRAP™
Absorbable Strap Fixation Device



ETHICON SECURESTRAP™ Open
Absorbable Strap Fixation Device

Codes	STRAP12, STRAP25	OPSTRAP20
Description	Single-use device with absorbable straps	Single-use open device with absorbable straps
Size	5mm	--
Strap Length	6.7mm	6.7mm
Strap Material	PLG Blend	PLG Blend
Strap Delivery Method	Expulsion	Expulsion
# of Tacks	12, 25	20
Qty./Box	6	6