






DEVICE SPECS

Ethicon Flat Hernia Mesh Portfolio

	 <p>ULTRAPRO ADVANCED™ Macroporous Partially Absorbable Mesh</p>	 <p>ULTRAPRO™ Macroporous Partially Absorbable Mesh</p>	 <p>PROLENE™ Soft Polypropylene Mesh</p>	 <p>PROLENE™ Polypropylene Mesh</p>	 <p>VICRYL™ (polyglactin 910) Woven & Knit Mesh</p>
Codes	UPA3612, UPA31515, UPA1530, UPA31015, UPA37615, UPA3030	UMM3, UMR3, UML1, UMS3	SPMXS, SPM3XL, SPMS, SPMXXL, SPMII, SPMLI, SPMH	PMII, PMXL, PMLK, PMSK, PMXS, PMH, PML	VKML, VKMM, VWML, VWMM
Description	Rectangular or square mesh	Rectangular or square mesh	Rectangular or square mesh	Rectangular, square, or large/ small pre-shaped keyhole mesh	Knitted or Woven Mesh
Size	6cm x 12cm, 15cm x 15cm, 15cm x 30cm, 10x15cm, 7.6cm x 15cm, 30cm x 30cm	15cm x 15cm, 7.6cm x 15cm, 30cm x 30cm, 6cm x 11cm	2.54cm x 10cm, 50cm x 50cm, 5cm x 10cm, 35.6cm x 30cm, 7.6cm x 15cm, 25cm x 25cm, 15cm x 15cm	8cm x 15cm, 5cm x 30cm, 6cm x13.7mm, 4.6cm x 10cm, 2.54cm x 10cm, 15cm x 15cm, 30cm x 30cm	30cm x 30cm, 15cm x 15cm
Material	Polypropylene, Poliglecaprone 25	Polypropylene, Poliglecaprone 25	Polypropylene	Polypropylene	Polyglactin 910
Qty./Box	1 (UPA1530, UPA3030) 3 (UPA3612, UPA31515, UPA31015, UPA37615)	1 (UML1) 3 (UMM3, UMR3, UMS3)	3 (SPM3XL, SPMXXL, SPMLI) 6 (SPMXS, SPMS, SPMII, SPMH)	3 (PML) 6 (all other codes)	3