


DEVICE SPECS

# Ethicon Hernia Mesh Devices



**ULTRAPRO™** Hernia System



**PROLENE™** Polypropylene Hernia System



**ULTRAPRO™** Plug

	<b>ULTRAPRO™</b> Hernia System	<b>PROLENE™</b> Polypropylene Hernia System	<b>ULTRAPRO™</b> Plug
<b>Codes</b>	UHSL6, UHSM6, UHSL, UHSM, UHSOV	PHSE6, PHSM6, PHSL6, PHSM, PHSE, PHSL	UPPM6, UPPS2, UPPS6, UPPM2, UPPL2, UPPL6
<b>Description</b>	Large, medium, or oval hernia system	Extended, medium, or large hernia system	Small, medium, large plugs
<b>Rim Size</b>	--	--	5cm
<b>Anchor Size</b>	--	--	3cm (UPPS2, UPPS6) 4cm (UPPM2, UPPM6) 5cm (UPPL2, UPPL6)
<b>Underlay Patch Size</b>	10cm (UHSL, UHSL6) 7.5cm (UHSM, UHSM6) 10cm x 12cm (UHSOV)	7.5cm (PHSM, PHSM6) 10cm (PHSL, PHSL6, PHSE, PHSE6)	--
<b>Onlay Patch Size</b>	6cm x 12cm	4.5cm x 10cm (PHSM, PHSM6, PHSL, PHSL6) 5.5cm x 12.8cm (PHSE, PHSE6)	--
<b>Connector Size</b>	--	1.9cm x 1.3cm	--
<b>Qty./Box</b>	3 (UHSL, UHSM, UHSOV) 6 (UHSL6, UHSM6)	3 (PHSM, PHSE, PHSL) 6 (PHSM6, PHSE6, PHSL6)	2 (UPPS2, UPPM2, UPPL2) 6 (UPPS6, UPPM6, UPPL6)