



DEVICE SPECS

LAPRA-TY™ Suture Clips

	 LAPRA-TY™ Absorbable Suture Clip Applier	 LAPRA-TY™ Absorbable Suture Clips
Codes	KA200	XC200
Description	Absorbable Suture Clip Applier	Absorbable Suture Clips Cartridge
Compatible Suture	--	2-0, 3-0, or 4-0 Coated VICRYL™ Suture
Length	33cm	--
Jaw Angle	22°	--
Shaft Rotation	360°	--
Clip Material	--	poly (p-dioxanone)
Reusable	Yes	No
Soft Tissue Approximation Period	--	Up to 14 days
Clips per Cartridge	--	6
Qty./Box	1	6