


DEVICE SPECS


Ethicon Linear Cutter Portfolio



ECHELON™ Linear Cutter



ETHICON™ Linear Cutter



PROXIMATE™ Linear Cutter

Codes	GLC60, GLC80, GLC100	NTLC55, NTLC75	TVC55*, TLC55, TCT55, TLC75, TCD75, TCT75, TLC10, TCT10
Associated Reloads	GLC60B, GLC60G, GLC80B, GLC80G, GLC100B, GLC100G	SR55, SR75	TVR55*†, TCR55, TRT55, TCR75, TRD75, TRT75, TCR10, TRT10
Description	Linear cutter with 3D/GST Reload	Selectable Linear Cutter with 3D Staples	Linear Cutter with Reload
Device Technology	Reloadable with locked halves	Selectable Staple Height	Cam Closure Mechanism
Size	60mm, 80mm, 100mm	55mm, 75mm	55mm, 75mm, 100mm
Staple Rows	4	6	4
Closed Staple Height	1.5mm (blue reloads), 2.0mm (green reloads)	1.5mm, 1.8mm, 2.0mm	1.0mm†, 1.5mm, 1.8mm, 2.0mm
Max Firings	8	12	8
Staple Technology	3D & GST	3D	B-Form
Qty./Box	3	3	3

*Vascular reload TVR55 only compatible with vascular stapler TVC55

† Vascular reload TVR55 closes to 1.0mm