


DEVICE SPECS

PROXIMATE® PPH Circular Stapler Set

ACCESSORIES INCLUDED



	PROXIMATE™ PPH Circular Stapler	Suture Threader	Purse-String Suture Anoscope	Circular Anal Dilator
Codes	PPH03	PPH03	PPH03	PPH03
Description	Hemorrhoidal Circular Stapler	Aids in pulling the sutures through the housing	Designed to facilitate purse-string suture application above the dentate line	Aid insertion of the Hemorrhoidal Circular Stapler into the anal canal.
Clear	--	--	Yes	Yes
Anvil	Fixed	--	--	--
Lumen Diameter	33mm	--	33mm	33mm
Staple Rows	2	--	--	--
Open Staple Leg Length	4mm	--	--	--
Closed Staple Height	0.75mm-1.5mm	--	--	--
Staple Technology	B-Form	--	--	--
Qty./Box	3	3	3	3