



DEVICE SPECS

# Ethicon Tissue Separating Hernia Mesh Portfolio



**ETHICON PHYSIOMESH™**  
Open Flexible Composite Mesh Device



**PROCEED™**  
Surgical Mesh



**PROCEED™**  
Ventral Patch

	<b>ETHICON PHYSIOMESH™</b> Open Flexible Composite Mesh Device	<b>PROCEED™</b> Surgical Mesh	<b>PROCEED™</b> Ventral Patch
Codes	OPHY1525, OPHY2536, OPHY2030, OPHY1215, OPHY1010, OPHY1520	PCDN1, PCDD1, PCDL1, PCDW1, PCDG1, PCDJ1, PCDH1, PCDT1, PCDM1, PCDR1, PCDB1	PVPM, PVPS
Description	Hexagonal mesh	Oval, rectangular, or square mesh	Circular mesh patch
Size	15cm x 25cm, 25cm x 36cm, 20cm x 30cm, 12cm x 15cm, 10cm round, 15cm x 20cm	10cm x 15cm, 10cm x 20cm, 30.5cm x 30.5cm, 25cm x 35cm, 15cm x 20cm, 20cm x 30cm, 20cm x 25cm, 26cm x 34cm, 15cm x 15cm, 7.5cm x 15cm, 5cm x 10cm	6.4cm x 6.4cm, 4.3cm x 4.3cm
Material	Polypropylene, Polydioxanone	Polypropylene	Polypropylene
Antibacterial	No	Yes	No
Qty./Box	1	1	2