

DEVICE SPECS

Ethicon Uterine & Pelvic Portfolio



GYNECARE INTERCEED™
Absorbable Adhesion Barrier



ARTISYN™ Y-Shaped Mesh



GYNECARE GYNEMESH™ PS
Nonabsorbable PROLENE™ Soft Mesh

Codes	4350, 4350XL	ARTY, ARTY5L	GPSL, GPSXL3
Description	Absorbable adhesion barrier	Partially absorbable mesh	Nonabsorbable bi-directional mesh
Size	7.6cm x 10.2cm (4350) 12.7cm x 15.2cm (4350XL)	27cm x 5cm	10cm x 15cm (GPSL) 25cm x 25cm (GPSXL3)
Material	Oxidized regenerated cellulose	Poliglecaprone 25/Polypropylene	Polypropylene
Surgical Approach	Open Gynecology	Sacrocolposuspension/ sacrocolpopexy (laparotomy or laparoscopic approach)	Sacrocolposuspension/ sacrocolpopexy (laparotomy or laparoscopic approach)
Storage	Store at or below 30°C	Store below 30°C	Store below 25°C
Qty./Box	10	1 (ARTY) 5 (ARTY5L)	1 (GPSL) 3 (GPSXL3)